



## **Scil Vet ABC - 1 Year Extended Warranty**

A great option to ensure 100% worry-free coverage!

As well as full coverage on all parts and labour, the extended warranty also includes one preventative maintenance.

### ***INCLUDED:***

#### **\* Preventative Maintenance**

- Complete change of o'rings & tubings
- Voltage measurement checks & adjustments
- All mechanical checkpoints tested
- Pressure checks on vacuum, needle sensor, carriage sensor
- Hemoglobin photometer calibration
- Latex gain adjustments (WBC, Gran, Lymph, RBC, PLT)
- Extensive cleaning inside & outside the analyzer
- Calibration & Control

**\* Complete warranty for one year on all parts & labour**

**\* Unlimited telephone technical support**

**\* Loaner unit**

***Call us today 1-866-382-6937***