



Scil Vet ABC - Preventative Maintenance

Vet Novations highly recommends that 15 months after the purchase of your Vet ABC a preventative maintenance is performed on an annual basis.

With the costs of parts, shipping and technician time all included, the annual maintenance is a smart choice to ensure that your clinic always guarantees patients the most accurate results.

Over time, as with any mechanical product, parts do wear due to everyday use. As stated below, the annual maintenance will include the replacement of these parts, so that you have a fresh start each year.

INCLUDED:

*** Preventative Maintenance**

- Complete change of o'rings & tubings
- Voltage measurement checks & adjustments
- All mechanical checkpoints tested
- Pressure checks on vacuum, needle sensor, carriage sensor
- Hemoglobin photometer calibration
- Latex gain adjustments (WBC, Gran, Lymph, RBC, PLT)
- Extensive cleaning inside & outside the analyzer
- Calibration & Control

*** Free loaner unit**

*** All shipping costs**

Call us today 1-866-382-6937