

Preventative Maintenance Log Book for scil Vet abc

Index: VCC = Vet abc concentrated cleaning (weekly)
 VCAL = Vet abc calibration (every 6 month)
 VBC = Vet abc blood control (daily)
 VCP = Vet abc change reagent pack
Year:

User List	Initials
A	
B	
C	

Vet Novations

Toll Free 1-866-382-6937

Day	Jan.	Feb.	March	April	May	June	July	Aug.	Sep.	Oct.	Nov.	Dec.
1												
2												
3												
4												
5												
6												
7												
8												
9												
10												
11												
12												
13												
14												
15												
16												
17												
18												
19												
20												
21												
22												
23												
24												
25												
26												
27												
28												
29												
30												
31												